


**Product Name: Recombinant Human F8 Protein**  
**Catalog #: PEH8105**



## Summary

<b>Name</b>	FA8_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P00451
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	33 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	NSSHLPAASY RKKDSGVQES SHFLQGAKKN NLSLAITLE MTGDQREVGS LGTSATNSVT YKKVENTVLP KPDLPKTSKG VELLPKVHIY QKDLFPTETS NGSPGHLDLV EGSLQGTGEG AIKWNEANRP GKVPFLRVAT ESSAKTPSKL LDPLAWDNHY GTQIPKEEWK SQEKSPEKTA FKKKDTILSL NACESNHAIA AINEGQNKPE IEVTWAKQGR TERLCSQNPV VLKRHQREIT RTTLQSDQEE IDYDDTISVE MKKEDFDIYD EDENQSPRSF

## image



## Background

<b>Alternative Names</b>	Coagulation factor VIII (Antihemophilic factor) (AHF) (Procoagulant component) [Cleaved into: Factor VIIIa heavy chain, 200 kDa isoform; Factor VIIIa heavy chain, 92 kDa isoform; Factor VIII B chain; Factor VIIIa light chain]
--------------------------	---

**Product Name: Recombinant Human F8 Protein**  
**Catalog #: PEH8105**

---



**Note**

For Research Use Only , Not for Diagnostic Use