


Product Name: Recombinant Human CUL1 Protein
Catalog #: PEH8113



Summary

Name	CUL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q13616
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RNGETINTRL ISGVVQSYVE LGLNEDDAFA KGPTLTVYKE SFESQFLADT ERFYTRESTE FLQQNPVTEY MKKAEARLLE EQRRVQVYLH ESTQDELARK CEQVLIEKHL EIFHTEFQNL LDADKNEDLG RMYNLSRIQ DGLGELKLL ETHIHNQGLA AIEKCGEAL NDPKMYVQTV

image



Background

Alternative Names Cullin-1 (CUL-1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CUL1 Protein
Catalog #: PEH8113

