


Product Name: Recombinant Human COX7A2L Protein
Catalog #: PEH8122



Summary

Name	COX7R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O14548
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKNKVPELQK FFQKADGVPV YLKRGLPDQM LYRTTMALTV GGTIYCLIAL YMASQPKNK

image



Background

Alternative Names	Cytochrome c oxidase subunit 7A-related protein, mitochondrial (COX7a-related protein) (Cytochrome c oxidase subunit VIIa-related protein) (EB1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use