


Product Name: Recombinant Human CYTIP Protein
Catalog #: PEH8124



Summary

Name	CYTIP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60759
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKQDNETFGF EIQSYRPQNQ NACSSEMFTL ICKIQEDSPA HCAGLQAGDV LANINGVSTE GFTYKQVVDL IRSSGNLLTI ETLNGTMILK RTELEAKLQV LKQTLKQKWV EYRSLQLQEH RLLHGDAANC PSLENMDLDE LSLFGPLPGP GPALVDRNRL SSESSCKSWL SSMTMDSEDG YQTCVSEDSS RGAFSRQTST DDECIPKEG

image



SDS

Background

Alternative Names	Cytohesin-interacting protein (Cytohesin binder and regulator) (CYBR) (Cytohesin-associated scaffolding protein) (CASP) (Cytohesin-binding protein HE) (Cbp HE) (Pleckstrin homology Sec7 and coiled-coil domains-binding protein)
--------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Product Name: Recombinant Human CYTIP Protein
Catalog #: PEH8124



Note

For Research Use Only , Not for Diagnostic Use