


Product Name: Recombinant Human DEAF1 Protein
Catalog #: PEH8128



Summary

Name	DEAF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75398
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGLAEAAVA AAAAVAAAAA AAAGGEAEEP VLSRDEEDSEE DADSEAERET PRVTAVAVMA AEPGHMDMGA EALPGPDEAA AAAAFAEVTT VTVANVGAAA DNVFTTSVAN AASISGHVLS GRTALQIGDS LNTEKATLIV VHTDGSIVET

image



Background

Alternative Names	Deformed epidermal autoregulatory factor 1 homolog (Nuclear DEAF-1-related transcriptional regulator) (NUDR) (Suppressin) (Zinc finger MYND domain-containing protein 5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DEAF1 Protein
Catalog #: PEH8128

