


Product Name: Recombinant Human DPP9 Protein
Catalog #: PEH8131



Summary

Name	DPP9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q86TI2
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GYAVVVIDGR GSCQRGLRFE GALKNQMGQV EIEDQVEGLQ FVAEKYGFID LSRVAIHGWS YGGFSLMGL IHKPQVKVA IAGAPVTVWM AYDTGYTERY MDVPENNHQHG YEAGSVALHV EKLPNENRNL LILHGFLDEN VHHFHTNFLV SQLIRAGKPY QLQIYPNERH SIRCPESEGH

image



Background

Alternative Names	Dipeptidyl peptidase 9 (DP9) (EC 3.4.14.5) (Dipeptidyl peptidase IV-related protein 2) (DPRP-2) (Dipeptidyl peptidase IX) (DPP IX) (Dipeptidyl peptidase-like protein 9) (DPLP9)
--------------------------	--

Note

Product Name: Recombinant Human DPP9 Protein
Catalog #: PEH8131



For Research Use Only , Not for Diagnostic Use