

Product Name: Recombinant Human RPN1 Protein
Catalog #: PEH8136



Summary

Name	RPN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P04843
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NRYKQSRDIS TLNSGKKSLE TEHKALTSEI ALLQSRLKTE GSDLCDRVSE MQKLDAQVKE LVLKSAVEAE RLVAGKLKGD TYIENEKLIS

image



Background

Alternative Names	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit 1 (Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 67 kDa subunit) (Ribophorin I) (RPN-I) (Ribophorin-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use