


Product Name: Recombinant Human NEURL1 Protein
Catalog #: PEH8139



Summary

Name	NEUL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O76050
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLCDLNVPGA DGDEAAPAAG CPIPQNSLNS QHSRALPAQL DGDLRFHALR AGAHVRILDE QTVARVEHGR DERALVFTSR PVRVAETIFV KVTRSGGARP GALSFGVTTC DPGTLRPADL PFSPEALVDR KEFWAVCRVP GPLHSGDILG LVVNADGELH LSHNGAAAGM QLCVDASQPL WMLFGLHGTI TQIRILGSTI

image



Background

Alternative Names	E3 ubiquitin-protein ligase NEURL1 (EC 2.3.2.27) (Neutralized-like protein 1A) (h-neu) (h-neutralized 1) (RING finger protein 67) (RING-type E3 ubiquitin transferase NEURL1)
--------------------------	---

Note

Product Name: Recombinant Human NEURL1 Protein
Catalog #: PEH8139



For Research Use Only , Not for Diagnostic Use