


Product Name: Recombinant Human UFL1 Protein
Catalog #: PEH8140



Summary

Name	UFL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94874
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESQSSHTGKK KPEISFMFQD EIEDFLRKHI QDAPEEFISE LAEYLIKPLN KTYLEVVRVS FMSSTTSASG TGRKRTIKDL

image



Background

Alternative Names E3 UFM1-protein ligase 1 (EC 2.3.2.-) (E3 UFM1-protein transferase 1) (Multiple alpha-helix protein located at ER) (Novel LZAP-binding protein) (Regulator of C53/LZAP and DDRGK1)

Note

For Research Use Only , Not for Diagnostic Use