


Product Name: Recombinant Human ETFA Protein
Catalog #: PEH8142



Summary

Name	ETFa_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	P13804
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HDVYKGLLPE ELTPLILATQ KQFNETHICA GASAFGKNLL PRVAAKLEVA PISDIIAISK PDLFVRLTIYA GNALCTVKCD EKVVFVSVRG TSFDAAATSG GSASSEKASS TSPVEISEWL DQKLTKSDRP ELTGAKVVVS GGRGLKSGEN FKLLYDLADQ LHAAVGASRA AVDAGFVPND MQVGQTGKIV APELYIAVGI SGAIQHLAGM KDSKTIVAIN KDPEAPIFQV ADYGIVADLF KVVPEMTEIL

image



Background

Alternative Names	Electron transfer flavoprotein subunit alpha, mitochondrial (Alpha-ETF)
--------------------------	---

Note

Product Name: Recombinant Human ETFA Protein
Catalog #: PEH8142



For Research Use Only , Not for Diagnostic Use