


Product Name: Recombinant Human GCKR Protein
Catalog #: PEH8154



Summary

Name	GCKR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14397
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IETPEPGKWE LSGYEAAPVI TEKSNPLTQD LDKADAENIV RLLGQCDAEI FQEEGQALST YQRLYSESL TTMVQVAGKV QEVLKEPDGG LVLSSGGGTS GRMAFLMSVS FNQLMKGLGQ KPLYTYLIAG GDRSVVASRE GTEDSALHGI EELKKVAAGK KRVIVIGISV GLSAPFVAGQ MDCCMNNTAV FLPVLVGFNP VSMARNDPIE DWSSTFRQVA ERMQKMQEKQ KAFVLNPAIG PEGLSGSSRM KGGSATKILL ETLLAAHKT VDQGIAASQR CLLEILRTFE RAHQVTYSQS

image



Background

Alternative Names Glucokinase regulatory protein (GKRP) (Glucokinase regulator)

Product Name: Recombinant Human GCKR Protein
Catalog #: PEH8154



Note

For Research Use Only , Not for Diagnostic Use