


Product Name: Recombinant Human MCF2L Protein
Catalog #: PEH8163



Summary

Name	MCF2L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15068
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LMQGSFSVWT DHKRGHTKVK ELARFKPMQR HFLHEKAVL FCKKRENGE GYEKAPSYSY KQSLNMAAVG ITENVKGDAL KFEIWYNARE EVYIVQAPTP EIKAAWVNEI RKVLTSQLQA CREASQHRAL EQSQLPLPA PTSTSPSRGN SRNIKKLEER KTDPLSLEGY VSSAPLTKPP EKGKGWSKTS HSLEAPEDDG GWSSAEEQIN SSDAEEDGGL GPKKLVPVKY TVVADHEKGG PDALRVRSGL VVELVQEGDE GLWYVRDPTT GKEGWVPASS

image



Background

Alternative Names	Guanine nucleotide exchange factor DBS (DBL's big sister) (MCF2-transforming sequence-like protein)
--------------------------	---

Product Name: Recombinant Human MCF2L Protein
Catalog #: PEH8163



Note

For Research Use Only , Not for Diagnostic Use