



Product Name: Recombinant Human GNG5 Protein
Catalog #: PEH8165



Summary

Name	GBG5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P63218
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSGSSVAAM KKVVQLRLE AGLNRVKVSQ AAADLKQFCL QNAQHDPLL GVSSSTNPFR PQQVC

image



Background

Alternative Names Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5

Note

For Research Use Only , Not for Diagnostic Use