


Product Name: Recombinant Human HDHD3 Protein
Catalog #: PEH8169



Summary

Name	HDHD3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q9BSH5
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAHRLQIRLL TWDVKDTLLR LRHPLGEAYA TKARAHGLEV EPSALEQGFR QAYRAQSHSF PNYGLSHGLT SRQWWLDVVL QTFHLAGVQD AQAVAPIAEQ LYKDFSHPCT WQVLDGAEDT LRECRTRGLR LAVISNFDRR LEGILGGLGL REHFDFVLTS EAAGWPKPDP RIFQEALRLA HMEPVVAHV GDNYLCDYQG PRAVGMHSFL VVGPQALDPV VRDSVPKEHI LPSLAHLLPA LDCLEGSTPG L

image



Background

Alternative Names Haloacid dehalogenase-like hydrolase domain-containing protein 3

Note

Product Name: Recombinant Human HDHD3 Protein
Catalog #: PEH8169



For Research Use Only , Not for Diagnostic Use