


Product Name: Recombinant Human HSPA13 Protein
Catalog #: PEH8171



Summary

Name	HSP13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P48723
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PHSSDTELPK DKLSSADDHR VNSGFGRGLS DKKSGESQVL FETEISRKLF DTLNEDLFQK ILVPIQQVLK EGHLEKTEID EVVLVGGSTR IPRIQVIQE FFGKDPNTSV DPDLAVVTGV

image



Background

Alternative Names	Heat shock 70 kDa protein 13 (Microsomal stress-70 protein ATPase core) (Stress-70 protein chaperone microsome-associated 60 kDa protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use