



Product Name: Recombinant Human HNF1A Protein
Catalog #: PEH8173



Summary

Name	HNF1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20823
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGGRGELAEI PNLGETRGS EDETDDGDED FTTPILKELE NLSPEEAHQ KAVVETLLQE DPWRVAKMVK SYLQQHNIPQ REVVDTTGLN QSHLSQHLNK GTPMKTQKRA ALYTWYVRKQ REVAQQFTHA GQGGLIEPT GDELPTKKGR RNRFKWGPAS QQILFQAYER QKNPSKEERE TLVEECNRAE CIQRGVSPSQ AQGLGSNLVT EVRVYNWFAN RRKEEAFRHK LAMDTYSGPP PGPGPGPALP AHSSPGLPPP

image



Background

Alternative Names	Hepatocyte nuclear factor 1-alpha (HNF-1-alpha) (HNF-1A) (Liver-specific transcription factor LF-B1) (LFB1) (Transcription factor 1) (TCF-1)
--------------------------	--

Product Name: Recombinant Human HNF1A Protein
Catalog #: PEH8173



Note

For Research Use Only , Not for Diagnostic Use