


Product Name: Recombinant Human ADHFE1 Protein
Catalog #: PEH8179



Summary

Name	HOT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8IWW8
Host	E.coli
Species	Human
Molecular Mass	65 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTDDYAFEMA VSNIRYGA AV TKEVGMDLKN MGAKNVCLMT DKNLSKLPPV QVAMDSL VKN GIPFTVYDNV RVEPTDSSFM EAIEFAQKGA FDAYVAVGGG STMDTCKAAN LYASSPHSDF LDYVSAPIGK GKPVSVPLKP LIAVPTTSGT GSETTGVAIF DYEHLKV KIG ITSRAIKPTL GLIDPLHTLH MPARVVANS G FDVLCHALE S

image



Background

Alternative Names	Hydroxyacid-oxoacid transhydrogenase, mitochondrial (HOT) (EC 1.1.99.24) (Alcohol dehydrogenase iron-containing protein 1) (ADHFe1) (Fe-containing alcohol dehydrogenase)
--------------------------	---

Product Name: Recombinant Human ADHFE1 Protein
Catalog #: PEH8179



Note

For Research Use Only , Not for Diagnostic Use