


Product Name: Recombinant Human KPNA6 Protein
Catalog #: PEH8180



Summary

Name	IMA7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60684
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VILNCSALPC LLHLLSPKE SIRKEACWTI SNITAGNRAQ IQAVIDANIF PVLIEILQKA EFRTRKEAAW AITNATSGGT PEQIRYLVSL GCIKPLCDLL TVMSKIVQV ALNGLLENILR LGEQEGKRSG SGVNPYCGLI EEAYGLDKIE FLQSHENQEI YQKAFDLIEH YFGVEDDDSS LAPQVDETQQ QFIFQQPEAP

image



Background

Alternative Names Importin subunit alpha-7 (Karyopherin subunit alpha-6)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KPNA6 Protein
Catalog #: PEH8180

