


**Product Name: Recombinant Human PPA2 Protein**  
**Catalog #: PEH8181**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | IPYR2_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | Q9H2U2   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 32 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.   |
| <b>Sequence</b>              | SPFHDIPLKV NSKEENGIPM KKARNDEYEN LFNMIVEIPR WTNAKMEIAT<br>KEPMNPIKQY VKDGKLRIVA NIFPYKGYIW NYGTLPQTWE DPHEKDKSTN<br>CFGDNDPIDV CEIGSKILSC GEVIHVKILG ILALIDEGET DWKLIAINAN<br>DPEASKFHDI DDVKKFKPGY LEATLNWFRL YKVPDGKPEN QFAFNGEFKN<br>KAFALEVIKS THQCWKALLM KKCNGGAINC TNVQISDSPF RCTQEEARSL<br>VESVSSSPNK |

## image



## Background

**Alternative Names** Inorganic pyrophosphatase 2, mitochondrial (EC 3.6.1.1) (Pyrophosphatase SID6-306) (Pyrophosphate phospho-hydrolase 2) (PPase 2)

**Product Name: Recombinant Human PPA2 Protein**  
**Catalog #: PEH8181**

---



**Note**

For Research Use Only , Not for Diagnostic Use