



**Product Name: Recombinant Human ILF3 Protein**  
**Catalog #: PEH8182**



## Summary

<b>Name</b>	ILF3_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q12906
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	25 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	KSPSKKKKKI QKKEEKAEPP QAMNALMRLN QLKPGQLQYKL VSQTGPVHAP IFTMSVEVDG NSFEEASGPSK KTAKLHVAVK VLQDMGLPTG AEGRDSSKGE DSAEETEAKP AVVAPAPVVE AVSTPSAAFP SDATEAQGPI LTKHGKNPVM ELNEKRRGLK YELISETGGS HDKRFVMEVE VDGQKFQGAG SNKKVAKAYA ALAALEKLFP DTPLALDANK

## image



## Background

<b>Alternative Names</b>	Interleukin enhancer-binding factor 3 (Double-stranded RNA-binding protein 76) (DRBP76) (M-phase phosphoprotein 4) (MPP4) (Nuclear factor associated with dsRNA) (NFAR) (Nuclear factor of activated T-cells 90 kDa) (NF-AT-90) (Translational control protein 80) (TCP80)
--------------------------	--

**Product Name: Recombinant Human ILF3 Protein**  
**Catalog #: PEH8182**

---



**Note**

For Research Use Only , Not for Diagnostic Use