


Product Name: Recombinant Human IARS2 Protein
Catalog #: PEH8186



Summary

Name	SYIM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NSE4
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NPEHVSRSIY VKFPLLKPSP KLASLIDGSS PVSILVWTTQ PWTIPANEAV CYMPESKYAV VKCSKSGDLY VLAADKVASV ASTLETTFET ISTLSGVDLE NGTCSHPLIP DKASPLL PAN HVTMAKGTGL VHTAPAHGME DYGVASQHNL PMDCLVDEDG VFTDVAGPEL

image



Background

Alternative Names	Isoleucine--tRNA ligase, mitochondrial (EC 6.1.1.5) (Isoleucyl-tRNA synthetase) (IleRS)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IARS2 Protein
Catalog #: PEH8186

