

Product Name: Recombinant Human KIF1C Protein
Catalog #: PEH8188



Summary

Name	KIF1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43896
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TESQIGPEEA MERLQETEKI IAEINETWEE KLRKTEALRM EREALLAEMG VAVREDGGTV GVFSPPKTPH LVNLTNEDPLM SECLLYHIKD GVTRVGGQVDM DIKLTGQFIR EQHCLFRSIP QPDGEVVVTL EPCEGAETVY NGKLVTEPLV LKSGNRIVMG KNHVFRFNHP EQARLERERG VPPPPGPPSE PVDWNFAQKE LLEQQGIDIK LEMEKRLQDL ENQYRKEKEE ADLLEQQRL YADSDSGDSS

image



Background

Alternative Names Kinesin-like protein KIF1C

Note

Product Name: Recombinant Human KIF1C Protein
Catalog #: PEH8188



For Research Use Only , Not for Diagnostic Use