


Product Name: Recombinant Human LRCH4 Protein
Catalog #: PEH8189



Summary

Name	LRCH4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75427
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPCPAEDLFP GHRYDGGLDS GFHSVDSGSK RWSGNESTDE FSELSFRISE LAREPRGPPE RKEDGSADGD PVQIDFIDSH VPGEDEERGT VEEQRPPELS PGAGDRERAP SSRREEPAGE ERRRPDTLQL WQERERRQQQ QSGAWGAPRK DSSLKPGLRA VVGGAAGVST QAMHNGSPKS SASQAGAAAG QGAPAPAPAS QEPLPIAGPA

image



Background

Alternative Names	Leucine-rich repeat and calponin homology domain-containing protein 4 (Leucine-rich repeat neuronal protein 4) (Leucine-rich neuronal protein)
--------------------------	--

Note

Product Name: Recombinant Human LRCH4 Protein
Catalog #: PEH8189



For Research Use Only , Not for Diagnostic Use