

Product Name: Recombinant Human SHOC2 Protein
Catalog #: PEH8190



Summary

Name	SHOC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UQ13
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKEKDPKVPS AKEREKEAKA SGGFGKESKE KEPKTKGKDA KDGKKDSSAA QPGVAFSVDN TIKRPNPAPG TRKKSSNAEV IKELNKREE NSMRLLSKR SIHILPSSIK ELTQLTELYL YSNKLQSLPA EVGCLVNLMT

image



Background

Alternative Names Leucine-rich repeat protein SHOC-2 (Protein soc-2 homolog) (Protein sur-8 homolog)

Note

For Research Use Only , Not for Diagnostic Use