


Product Name: Recombinant Human LONP1 Protein
Catalog #: PEH8191



Summary

Name	LONM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P36776
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RNPVFPRFIK IIEVKNNKLV ELLRRKVRLA QPYVGVFLKR DDSNESDVVE SLDEIYHTGT FAQIHEMQDL GDKLRMIVMG HRRVHISRQL EVEPEEPEAE NKKHKPRRKSK RKGKEAEDEL SARHPAELAM EPTPELPAEV LMVEVENVVH EDFQVTEEVK ALTAEIVKTI RDIIALNPLY RESVLQMMQA GQRVVDNPIY LSDMGAALTG AESHELQDVL EETNIPKRLY KALSLLKKEF ELSKLQQLRG

image



Background

Alternative Names	Lon protease homolog, mitochondrial (EC 3.4.21.53) (LONHs) (Lon protease-like protein) (LONP) (Mitochondrial ATP-dependent protease Lon) (Serine protease 15)
--------------------------	---

Product Name: Recombinant Human LONP1 Protein
Catalog #: PEH8191



Note

For Research Use Only , Not for Diagnostic Use