



Product Name: Recombinant Human MVP Protein
Catalog #: PEH8195



Summary

Name	MVP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14764
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTVQDTEAHV PDVHEEVLGV VPITTLGPHN YCVILDPVGP DGKNQLGQKR VVKGEKSFFL QPGEQLEQGI QDVYVLSEQQ GLLLRALQPL EEGEDEEKVS HQAGDHWLIR GPLEYVPSAK VEVVEERQAI PLDENEGIYV QDVKTGKVRA VIGSTYMLTQ DEVLWEKELP PGVEELLNKG

image



Background

Alternative Names Major vault protein (MVP) (Lung resistance-related protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MVP Protein
Catalog #: PEH8195

