


Product Name: Recombinant Human MARK3 Protein
Catalog #: PEH8199



Summary

Name	MARK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P27448
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PHIGNYRLK TIGKGNFAKV KLARHILTGR EVAIKIIDKT QLNPTSLQKL FREVRIMKIL NHPNIVKLFE VIETEKTLYL IMEYASGGEV FDYLVAHGRM KEKEARSKFR QIVSAVQYCH QKRIVHRDLK AENLLLDADM NIKIADFGFS NEFTVGGKLD TFCGSPPYAA PELFQGKKYD GPEVDVWSLG VILYTLVSGS LPFDGQNLKE LRERVLRGKY RIPFYMSTDC ENLLKRFLVL NPIKRGTLQ IMKDRWINAG HEEDELKPFV EPELDSQDK RIDIMVGMGY SQEEIQESLS KMKYDEITAT YLLGRKSSE LDASDSSSSS

image



Background

Alternative Names	MAP/microtubule affinity-regulating kinase 3 (EC 2.7.11.1) (C-TAK1) (cTAK1) (Cdc25C-associated protein kinase 1) (ELKL motif kinase 2) (EMK-2) (Protein
--------------------------	---

Product Name: Recombinant Human MARK3 Protein
Catalog #: PEH8199



kinase STK10) (Ser/Thr protein kinase PAR-1) (Par-1a) (Serine/threonine-protein kinase p78)

Note

For Research Use Only , Not for Diagnostic Use