


Product Name: Recombinant Human TIMM50 Protein
Catalog #: PEH8208



Summary

Name	TIM50_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3ZCQ8
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPCLLPDPLQ EPYYQPPYTL VLELTGVLLH PEWSLATGWR FKKRPGIETL FQQLAPLYEI VIFTSETGMT AFPLIDSVDP HGFISYRLFR DATRYMDGHH VKDISCLNRD PARVVVVDCCK KEAFRLQPYN GVALRPWDGN SDDRVLDDL AFLKTIALNG VEDVRTVLEH YALEDDPLAA

image



Background

Alternative Names Mitochondrial import inner membrane translocase subunit TIM50

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TIMM50 Protein
Catalog #: PEH8208

