


Product Name: Recombinant Human MAP3K9 Protein
Catalog #: PEH8212



Summary

Name	M3K9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P80192
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PLSPCTHNPL VNVVERFKR DPNQSLTPTH VTLTTPSQPS SHRRTPSDGA LKPETLLASR SPSSNGLSPS PGAGMLKTPS PSRDGGEFPR LPDPNVVFPP TPRRWNTQQD STLERPKTLE FLPRPRPSAN RQRLDPWWFV SPSHARSTSP

image



Background

Alternative Names	Mitogen-activated protein kinase kinase kinase 9 (EC 2.7.11.25) (Mixed lineage kinase 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use