


## Summary

<b>Name</b>	MRP6_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O95255
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	34 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GSAAGKDCIT IHSATFAWSQ ESPCLHRIN LTVPQGCLLA VVGPGVAGKS SLLSALLGEL SKVEGFVSIE GAVAYVPQEA WVQNTSVVEN VCFGQELDPP WLERVLEACA LQPDVDSFPE GIHTSIGEQG MNLSGGQKQR LSLARAVYRK AAVYLLDDPL AALDAHVGQH VFNQVIGPGG LLQGTTRILV THALHILPQA DWIIVLANGA IAEMGSYQEL LQRKGALMCL LDQARQPGDR GEGETEPGTS TKDPRGTSAG RPELRRERS IKSVPKDR TSEAQTEVPL DDPDRAGWPA

## image



## Background

<b>Alternative Names</b>	Multidrug resistance-associated protein 6 (ATP-binding cassette sub-family C member 6) (Anthracycline resistance-associated protein) (Multi-specific organic anion transporter E) (MOAT-E)
--------------------------	--

**Product Name: Recombinant Human ABCC6 Protein**  
**Catalog #: PEH8214**

---



**Note**

For Research Use Only , Not for Diagnostic Use