


Product Name: Recombinant Human ME2 Protein
Catalog #: PEH8217



Summary

Name	MAOM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P23368
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGTAAVALAG LLAAQKVISK PISEHKILFL GAGEAALGIA NLIVMSMVEN GLSEQEAQKK IWMFDKYGLL VKGRKAKIDS YQEPFTHSAP ESIPDTFEDA VNILKPSTII GVAGAGRLFT PDVIRAMASI NERPVIFALS NPTAQAECTA EEAYTLTEGR CLFASGSPFG PVKLT DGRVF TPGQGNNVYI FPGVALAVIL CNTRHISDSV FLEAAKALTS QLTDEELAQG RLYPPLANIQ EVSINIAIKV TELYANKMA FRYPEPEDKA KYVKERTWRS

image



Background

Alternative Names	NAD-dependent malic enzyme, mitochondrial (NAD-ME) (EC 1.1.1.38) (Malic enzyme 2)
--------------------------	---

Product Name: Recombinant Human ME2 Protein
Catalog #: PEH8217



Note

For Research Use Only , Not for Diagnostic Use