


Product Name: Recombinant Human NASP Protein
Catalog #: PEH8221



Summary

Name	NASP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49321
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQETKEAQLY AAQAHLKLGE VSVESYNYVQ AVEEFQSCLN LQEQYLEAHD RLLAETHYQL GLAYGYNSQY DEAVAQFSKS IEVIENRMAV LNEQVKEAEG SSAEYKKEIE ELKELLPEIR EKIEDAKESQ RSGNVLAELAL KATLVESSTS GFTPGGGGSS VSMIASRKPT DGASSSNCVT DISHLVRKKR KPEEESPRKD DAKKAKQEPE VNGGSGDAVP SGNEVSENME EEAENQAESR

image



Background

Alternative Names Nuclear autoantigenic sperm protein (NASP)

Note

Product Name: Recombinant Human NASP Protein
Catalog #: PEH8221



For Research Use Only , Not for Diagnostic Use