


**Product Name: Recombinant Human NFATC4 Protein**  
**Catalog #: PEH8222**



## Summary

<b>Name</b>	NFAC4_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q14934
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	40 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SVAPPGSRK EVAGMDYLAV PSPLAWSKAR IGGHSPIFRT SALPPLDWPL PSQYEQLER IEVQPRHHR AHYETEGSRG AVKAAPGGHP VVKLLGYSEK PLTLQMFIGN ADERNLRPHA FYQVHRITGK MVATASYEAV VSGTKVLEMT LLPENNMAAN IDCAGILKLR NSDIELRKGE TDIGRKNTRV RLVFRVHVVPQ GGGKVVSVQA ASVPIECSQR SAQELPQVEA YSPSACSVRG GEELVLTGSN FLPDSKVVFI ERGPDGKLQW EEEATVNRLQ SNEVTLTLTV PEYSNKRVS RVQVYFYVSN GRRKRSPTQS FRFLPVICKE EPLPDSSLRG FPSASATPFG

## image



SDS

## Background

**Alternative Names** Nuclear factor of activated T-cells, cytoplasmic 4 (NF-ATc4) (NFATc4) (T-cell transcription factor NFAT3) (NF-AT3)

**Product Name: Recombinant Human NFATC4 Protein**  
**Catalog #: PEH8222**

---



**Note**

For Research Use Only , Not for Diagnostic Use