



Product Name: Recombinant Human NUCB1 Protein
Catalog #: PEH8226



Summary

Name	NUCB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q02818
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPLERGAPNK EETPATESPD TGLYYHRYLQ EVIDVLETDG HFREKLQAAN AEDIKSGKLS RELDFVSHHV RTKLDELKRQ EVSRLRMLLK AKMDAEQDPN VQVDHLNLLK QFEHLDPQNN HTFEARDLEL LIQTATRDLA QYDAAHHEEF KRYEMLKEHE RRRYLESIGE EQRKEAERKL EEQRRRHREH PKVNVPGSQA QLKEVWEELD GLDPNRFNPK TFFILHDINS DGVLDEQELE ALFTKELEKV YDPKNEEDDM REMEEERLRM REHVMKNVDT NQDRLVTLEE FLASTQRKEF GDTGEGWETV EMHPAYTEEE LRRFEEELAA REAELNAKAQ RLSQETEALG RSQGRLEAQK RELQQAVLHM EQRKQQQQQQ QGHKAPAAHP EGQLKFHPDT DDVPVPAPAG DQKEVDTSEK KLLERLPEVE VPQHL

image



Background

Product Name: Recombinant Human NUCB1 Protein
Catalog #: PEH8226



Alternative Names

Nucleobindin-1 (CALNUC)

Note

For Research Use Only , Not for Diagnostic Use