


Product Name: Recombinant Human PEX5 Protein
Catalog #: PEH8231



Summary

Name	PEX5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P50542
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGTSDAWVDQ FTRPVNTSAL DMEFERAKSA IESDVDFWDK LQAELEEMAK RDAEAHPWLS DYDDLTSATY DKG YQFEEEN PLRDHPQPFE EGLRRRLQEGD LPNAVLLFEA AVQQDPKHME AWQYL GTTQA ENEQELLAIS ALRRCLELKP DNQTALMALA

image



Background

Alternative Names	Peroxisomal targeting signal 1 receptor (PTS1 receptor) (PTS1R) (PTS1-BP) (Peroxin-5) (Peroxisomal C-terminal targeting signal import receptor) (Peroxisome receptor 1)
--------------------------	---

Note

Product Name: Recombinant Human PEX5 Protein
Catalog #: PEH8231



For Research Use Only , Not for Diagnostic Use