


Product Name: Recombinant Human PSD Protein
Catalog #: PEH8232



Summary

Name	PSD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A5PKW4
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPSAATRLSQ EEQVRTHEAK LKAMASELRE HRAAQLGKKG RGKEAEEQRQ KEAYLEFEKS RYSTYAALLR VKLKAGSEEL DAVEAALAQA GSTEDGLPPS HSSPSLQPKP SSQPRAQRHS SEPRPGAGSG

image



Background

Alternative Names	PH and SEC7 domain-containing protein 1 (Exchange factor for ADP-ribosylation factor guanine nucleotide factor 6) (Exchange factor for ARF6) (Exchange factor for ARF6 A) (Pleckstrin homology and SEC7 domain-containing protein 1)
--------------------------	--

Note

Product Name: Recombinant Human PSD Protein
Catalog #: PEH8232



For Research Use Only , Not for Diagnostic Use