


Product Name: Recombinant Human PIK3C2A Protein
Catalog #: PEH8233



Summary

Name	P3C2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00443
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KECPSSHPEP TRAKDVKKEE ALQMEAEALA KLQKDRQVTD NQRGFELSSS TRKKAQVYNK QDYDLMVFPE SDSQKRALDI DVEKLTQAEL EKLLDDSF TKKTPVLPVT PILSPSFAQ

image



Background

Alternative Names	Phosphatidylinositol 4-phosphate 3-kinase C2 domain-containing subunit alpha (PI3K-C2-alpha) (PtdIns-3-kinase C2 subunit alpha) (EC 2.7.1.154) (Phosphoinositide 3-kinase-C2-alpha)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PIK3C2A Protein
Catalog #: PEH8233

