


Product Name: Recombinant Human PIK3C2G Protein
Catalog #: PEH8234



Summary

Name	P3C2G_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75747
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PNESHEKQYE HQEFLFVNQP HSSSQVSLGF DQIVDEISGK IPHYESEIDE NTFFVPTAPK WDSTGHSLNE AHQISLNEFT

image



Background

Alternative Names	Phosphatidylinositol 4-phosphate 3-kinase C2 domain-containing subunit gamma (PI3K-C2-gamma) (PtdIns-3-kinase C2 subunit gamma) (EC 2.7.1.154) (Phosphoinositide 3-kinase-C2-gamma)
--------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

For Research Use Only , Not for Diagnostic Use