


Product Name: Recombinant Human PLS1 Protein
Catalog #: PEH8236



Summary

Name	PLSI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14651
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EELMKLSPEE LLLRWVNYHL TNAGWHTISN FSQDIKDSRA YFHLLNQIAP KGGEDGPAIA IDLSGINETN DLKRAGLMLQ EADKLGCKQF VTPADVVSNG PKLNLAFFAN LFNTYPLCHK PNNNDIDMNL LEGESKEERT FRNWMNSLGV NPYINHLYSD LADALVIFQL YEMIRVPVNW SHVKNPPYPA LGGNMCKIEN CNYAVELGKN KAKFSLVGIA GQDLNEGNST LTLALVWQLM

image



Background

Alternative Names Plastin-1 (Intestine-specific plastin) (I-plastin)

Note

Product Name: Recombinant Human PLS1 Protein
Catalog #: PEH8236



For Research Use Only , Not for Diagnostic Use