


Product Name: Recombinant Human KCNK10 Protein
Catalog #: PEH8238



Summary

Name	KCNKA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P57789
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIHDKLQRAA TIRSMERRRL GLDQRAHSLD MLSPEKRSVF AALDTGRFKA SSQESINNRPNLRLKGPEQLNKHGQGASEDNIINKFGSTSRLTKRKNKD

image



Background

Alternative Names	Potassium channel subfamily K member 10 (Outward rectifying potassium channel protein TREK-2) (TREK-2 K(+)) channel subunit)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use