


Product Name: Recombinant Human FAM114A1 Protein
Catalog #: PEH8257



Summary

Name	NXP20_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IWE2
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEFIGKTMN VLAESDPGFK RTKTMERTV SLSQMLREAK EKEKQLAQQ LTMERTAHYG MLFDEYQGLS HLEALEILSN ESESKVQSFL ASLDGEKLEL LKNDLISIKD IFAAKELENE ENQEEQGLEE KGEEFARMLT ELLFELHVAA TPDKLNKAMK RAHDWVEEDQ TVVSVDVAKV SEETKKEEK EEKSQDPQED KKEEKTKTI EEVYMSSIES LAEVTARCIE QLHKVAELIL HGQEEKPAQ DQAKVLIKLT TAMCNEVASL SKKFTNSLTT VGSNKKAEVL NPMISSVLE GCNSTTYIQD AFQLLLPVLQ VSHIQTSCLK

image



Background

Alternative Names Protein NOXP20 (Nervous system overexpressed protein 20) (Protein FAM114A1)

Product Name: Recombinant Human FAM114A1 Protein
Catalog #: PEH8257



Note

For Research Use Only , Not for Diagnostic Use