


Product Name: Recombinant Human SEC24C Protein
Catalog #: PEH8262



Summary

Name	SC24C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53992
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPAFIFMIDV SYNAIRTGLV RLLCEELKSL LDFLPREGGA EESAIRVGFV TYNKLHFN VKSSLAQPQM MVVSDVADMF VPLLDGFLVN VNESRAVITS LLDQIPEMFA DTRETETVFV PVIQAGMEAL KAAECAGKLF LFHTSLPIAE APGKLKNRDD RKLINTDKEK TLFQPQTGAY QTLAKECVAQ GCCVDLFLFP NQYVDVATLS VVPQLTGGSV YKYASFQVEN DQERFLSDLR

image



Background

Alternative Names Protein transport protein Sec24C (SEC24-related protein C)

Note

Product Name: Recombinant Human SEC24C Protein
Catalog #: PEH8262



For Research Use Only , Not for Diagnostic Use