


Product Name: Recombinant Human PNPO Protein
Catalog #: PEH8264



Summary

Name	PNPO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NVS9
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	THLTSLDPVK QFAAWFEEAV QCPDIGEANA MCLATCTRDG KPSARMLLLK GFGKDGFRFF TNFESRKGKE LDSNPFASLV FYWEPLNRQV RVEGPVKKLP EEEAE CYFHS RPKSSQIGAV VSHQSSVIPD REYLRKKNEE LEQLYQDQEV PKPKSWG GYV LYPQVMEFWQ GQTNRLHDRI VFRRGLPTGD SPLGPMTHRG

image



Background

Alternative Names Pyridoxine-5'-phosphate oxidase (EC 1.4.3.5) (Pyridoxamine-phosphate oxidase)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PNPO Protein
Catalog #: PEH8264

