



Product Name: Recombinant Human RAB3GAP1 Protein
Catalog #: PEH8266



Summary

Name	RB3GP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15042
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKETDKEKGE VGKSWDSWSD SEEEFFECLS DTEELKGNGQ ESGKKGGPKE MANLRPEGRL YQH GKLTLLH NGEPLYIPVT QEPAPMTEDL LEEQSEVLAK LG TSAEGAHL RARMQSACLL SDMESFKAAN PGCSLEDFVR WYSPRDYIEE EVIDEKGNNV LKGELSARMK IPSNMWVEAW ETAKPIPARR QRRLFDDTRE

image



Background

Alternative Names	Rab3 GTPase-activating protein catalytic subunit (RAB3 GTPase-activating protein 130 kDa subunit) (Rab3-GAP p130) (Rab3-GAP)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RAB3GAP1 Protein
Catalog #: PEH8266

