


Product Name: Recombinant Human RGS19 Protein
Catalog #: PEH8271



Summary

Name	RGS19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49795
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPTPHEAEKQ ITGPEEADRP PSMSSHD TAS PAAPSRNPCC LCWCCCCSCS WNQERRRAWQ ASRESKLQPL PSCEVCATPS PEEVQSWAQS FDKLMHSPAG RSVFRAFLRT EYSEENMLFW LACEELKAEA NQHVVDEKAR LIYEDYVSIL SPKEVSLDSR VREGINKKMQ EPSAHTFDDA QLQIYTLMHR DSYPRFLSSP TYRALLQGP SQSSSEA

image



Background

Alternative Names	Regulator of G-protein signaling 19 (RGS19) (G-alpha-interacting protein) (GAIP)
--------------------------	--

Note

Product Name: Recombinant Human RGS19 Protein
Catalog #: PEH8271



For Research Use Only , Not for Diagnostic Use