

Product Name: Recombinant Human RLBP1 Protein
Catalog #: PEH8273



Summary

Name	RLBP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12271
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLEQLTTKDH GPVFGPMSQL PRHTLQKAKD ELNEREETRE EAVRELQEMV QAQAASGEEL AVAVAERVQE KDSGFFLRFI RARKFNVGRA YELLRGYVNF RLQYPELFDS LSPEAVRCTI EAGYPGVLS RDYGRVVML FNINWQSQE ITFDEILQAY CFILEKLEN EETQINGFCI IENFKGFTMQ QAASLRSD RKMVDMLQDS FPARFKAIHF IHQPWYFTTT YNVVKPFLKS KLLERVFVHG DDLSGFYQEI DENILPSDFG

image



Background

Alternative Names Retinaldehyde-binding protein 1 (Cellular retinaldehyde-binding protein)

Product Name: Recombinant Human RLBP1 Protein
Catalog #: PEH8273



Note

For Research Use Only , Not for Diagnostic Use