


Product Name: Recombinant Human RBL2 Protein
Catalog #: PEH8275



Summary

Name	RBL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08999
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKGKEENLTG FLEPGNFGES FKAINKAYEE YVLSVGNLDE RIFLGEDAE EIGTLRCLN AGSGTETAER VQMKNILQQH

image



Background

Alternative Names	Retinoblastoma-like protein 2 (130 kDa retinoblastoma-associated protein) (p130) (Retinoblastoma-related protein 2) (RBR-2) (pRb2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use