


Product Name: Recombinant Human ARHGEF7 Protein
Catalog #: PEH8276



Summary

Name	ARHG7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14155
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GFDTTAINKS YYNVVLQNIL ETENEYSKEL QTVLSTYLRP LQTSEKLSSA NISYLMGNLE EICSFQQMLV QSLEECTKLP EAQQRVGGCF LNLMPQMKTL YLTTCANHPS AVNVLTEHSE ELGEFMETKG ASSPGILVLT TGLSKPFMRL DKYPTLLKEL ERHMEDYHTD RQDIQKSMAA FKNLSAQCQE VRKRKELELQ

image



Background

Alternative Names	Rho guanine nucleotide exchange factor 7 (Beta-Pix) (COOL-1) (PAK-interacting exchange factor beta) (p85)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARHGEF7 Protein
Catalog #: PEH8276

