


**Product Name: Recombinant Human PPIL2 Protein**  
**Catalog #: PEH8279**



## Summary

<b>Name</b>	PPIL2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q13356
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	23 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	YVRLHTNKGD LNLELHCDLT PKTCENFIRL CKKHYYDGTI FHRSIRNFVI QGGDPTGTGT GGESYWGKPF KDEFRPNLSH TGRGILSMAN SGPNSNRSQF FITFRSCAYL DKKHTIFGRV VGGFDVLTAM ENVESDPKTD RPKEEIRIDA TTVFVDPYEE ADAQIAQERK TQLKVAPETK VKSSQPQAGS

## image



## Background

**Alternative Names** RING-type E3 ubiquitin-protein ligase PPIL2 (EC 2.3.2.27) (CYC4) (Cyclophilin-60) (Cyclophilin-like protein Cyp-60) (Cyp60) (hCyP-60) (Probable inactive peptidyl-prolyl cis-trans isomerase-like 2) (PPIase) (Rotamase PPIL2)

**Product Name: Recombinant Human PPIL2 Protein**  
**Catalog #: PEH8279**



---

**Note**

For Research Use Only , Not for Diagnostic Use