


Product Name: Recombinant Human RFOX2 Protein
Catalog #: PEH8280



Summary

Name	RFOX2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43251
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQEPTTTPDA MVQPFTTIPF PPPPQNGIPT EYGVPHQTQDY AGQTGEHNLT LYGSTQAHGE QSSNSPSTQN GSLTTEGGAQ TDGQQSQTQS SENSESKSTP KRLHVSNIPI RFRDPLRQM FGQFGKILDV EIIIFNERGSK GFGFVTFENS ADADRAREKL HGTVVVEGRKI EVNNATARVM TNKKMVTPIYA NGWKLSPVVG AVYGPELYAA SSFQADVSLG NDAAVPLSGR GGINTYIPLI SLPLVPGFPY PTAATTAAAF RGAHLRGRGR TVYGAVRAVP PTAIPAYPGV VYQDGFYCAD

image



Background

Alternative Names	RNA binding protein fox-1 homolog 2 (Fox-1 homolog B) (Hexaribonucleotide-binding protein 2) (RNA-binding motif protein 9) (RNA-binding protein 9) (Repressor of tamoxifen transcriptional activity)
--------------------------	--

Product Name: Recombinant Human RBFOX2 Protein
Catalog #: PEH8280



Note

For Research Use Only , Not for Diagnostic Use